

MagCore® Plasma DNA Extraction Kit (1.2 ml)

For extraction of free circulating DNA from human plasma or serum.

Applicable Models : HF16, Compact, HF48, Super, HF16 Plus, Plus II

Cartridge Code 105

Cat.No. MPD1200

Kit Contents


Check that the following parts are included in addition to the main unit:

Cat.No. MPD1200 Contents:	
Pre-filled Cartridge Reagent.....	96 pcs.
Pipet Tip plus Holder Set.....	100 sets.
Sample Tube.....	100 pcs.
Elution Tube.....	100 pcs.
Proteinase K(11 mg).....	2 pcs.
PK Storage Buffer.....	2 pcs.

Storage and Stability:

1. This kit should be stored at room temperature.
2. Proteinase K should be stored at 2-8 °C upon arrival.
3. Shelf Life : 18 Months.

Cartridge Contents :



Description

MagCore® Plasma DNA Extraction Kit is designed for purification of DNA from 1.2 ml of serum, plasma, cell-free body fluids by using MagCore® auto-extraction instrument. With all the kit components of plastic consumables are DNase/ RNase-Free pretreated and individual processing track for each loaded sample, this system eliminates all possible cross contamination between samples. Built-in protocol with flexibility in sample source volumes, plasma DNA can be extracted using this kit in a fast and economical way.

Applications

The purified total nucleic acid is suitable for highly sensitive and quantitative PCR. MagCore® Plasma DNA Extraction Kit has been proven with various genomic analyses as downstream applications.

Preparation Before Using

1. Add 1.1 ml PK Storage Buffer to the Proteinase K tube and mix by vortexing. Store prepared Proteinase K (10 mg/ml) at 2-8 °C

Protocol

1. Pipet 20 µl proteinase K (10 mg/ml) into the MagCore® Sample Tubes.
2. Add 1200 µl of serum, plasma, cell-free body fluids into the prepared Sample Tube.
3. After Proteinase K and Plasma mixing, stand for 10-20 min at room temperature, then centrifuge at 14,000 rpm for 5 min and transfer clear plasma to new tube.
4. Put the prepared Sample Tube into the correct well of T-Rack. (see page 3-10)
5. Put Elution Tube and Tip Plus Holder Set (HF16, Compact) / Pipette Tip (Super, Plus) into the correct wells of the T-Rack. (see page 3-10)
6. Run Code 105 program at MagCore®.

MagCore® Circulating DNA Large Volume Kit (4 ml)

For extraction of free circulating DNA from human plasma or serum.

Applicable Models : HF16, Compact, HF48, Super, HF16 Plus, Plus II

Cartridge Code 115

Cat.No. MPD4000-01 // MPD4000-03

Kit Contents

Check that the following parts are included in addition to the main unit:

Cat.No. MPD4000-01 Contents:

Pre-filled Cartridge Reagent.....	24 pcs.
Pipet Tip plus Holder Set.....	25 sets.
5ml Sample Tube.....	25 pcs.
Elution Tube.....	25 pcs.
Proteinase K(11mg).....	3 pcs.
PK Storage Buffer.....	3 pcs.


Cat.No. MPD4000-03 Contents:

Pre-filled Cartridge Reagent.....	96 pcs.
Pipet Tip plus Holder Set.....	100 sets.
5ml Sample Tube.....	100 pcs.
Elution Tube.....	100 pcs.
Proteinase K(11mg).....	10 pcs.
PK Storage Buffer.....	10 pcs.

Storage and Stability:

1. This kit should be stored at room temperature.
2. Proteinase K should be stored at 2-8 °C upon arrival.
3. Shelf Life: 18 Months.

Cartridge Contents :



Description

MagCore® Circulating DNA large volume kit is designed for purification of DNA from 3 ml or 4 ml of serum, plasma, cell-free body fluids by using MagCore® auto-extraction instrument. With all the kit components of plastic consumables are DNase/RNase-Free pretreated, and individual processing track for each loaded samples, this system eliminates all possible cross contamination between samples. Built-in protocol with flexibility in sample source volumes, plasma DNA can be extracted using this kit in a fast and economical way.

Applications

The purified total nucleic acid is suitable for highly sensitive and quantitative PCR. MagCore® Circulating DNA Large Volume Kit has been proven with various genomic analyses as downstream applications.

Preparation Before Using

Add 1.1 ml PK Storage Buffer to the Proteinase K tube and mix by vortexing. Store prepared Proteinase K (10 mg/ml) at 2-8 °C

Protocol

1. Add 4 ml of serum, plasma, cell-free body fluids (If the volume is less than 4 ml, add to 4 ml with 1X PBS.) into a DNase-free 15 ml tube (not provided).
2. Add 100 µl proteinase K (10 mg/ml) into 15 ml tube and mix by vortexing.
3. After Proteinase K and Plasma mixing, stand for 10-15 min at room temperature, then centrifuge at 14,000 rpm for 5 min and transfer clear plasma to 5 ml Sample Tube.
4. Put the prepared Sample Tube into the well 6 of the T-Rack.
5. Put the Pipette Tip into the well 3 of the T-Rack and the Elution Tube into the well 5 of the T-Rack.
6. Run Code 115 program at MagCore®.